

Mixing the Alpron solution to get a 1% concentration is crucial. Please refer to the chart below.

Mixed Solution(ml)	Alpron	Water	Mixed Solution (Ltr)
500 ml	5 ml	495 ml	0.5 Ltr
750 ml	7.5 ml	742.5 ml	0.75 Ltr
1000 ml	10 ml	990 ml	1.0 Ltr
1500 ml	15 ml	1485 ml	1.5 Ltr
2000 ml	20 ml	1980 ml	2.0 Ltr
2500 ml	25 ml	2475 ml	2.5 Ltr
3000 ml	30 ml	2970 ml	3.0 Ltr
3500 ml	35 ml	3465 ml	3.5 Lt
4000 ml	40 ml	3960 ml	4.0 Ltr
4500 ml	45 ml	4455 ml	4.5 Ltr
5000 ml	50 ml	4950 ml	5.0 Ltr